


Sealing Film

Sealing Film provides a highly effective seal, preventing evaporation and eliminating contamination between wells and are designed for applications ranging from PCR and Real-Time PCR to ELISA and cell culture.



- 2 types: standard, optical clear sealing film
- Pre-cut, plate-sized sheets
- Standard film is used for short-term storage and transportation of samples to ensure tightness and reduce evaporation of samples
- Optical clear sealing film (no auto-fluorescence, high transparency) allows for optical analysis during real-time PCR
- Pre-cut, plate-sized sheets
- Pressure sensitive film, non - adhesive to skin and gloves, convenient for experimental operation, does not affect optical analysis
- Not react with experimental samples ensure reliable results
- DNase-/RNase-free
- Nonsterile