

Centrifuge tubes

LABSELECT 15ml, 50ml centrifuge tubes are made of high-quality polypropylene (PP) and are mainly used for sample storage and centrifugation. These leak-proof conical-bottom centrifuge tubes prevent breakage and leakage during high-speed centrifugation. These disposable conical-bottom tubes feature black printed graduations and smooth inner walls for easy filling and sample preparation



Tight sealing screw cap

- Black printed graduations and a large white marking area
- Maximum RCF is 12,000xg
- Nonpyrogenic
- Noncytotoxic
- DNase-/RNase-free
- Sterile