



High Speed Centrifuge

Model: L-450C

Brush Less Motor

Max Speed 16000r/min

Max RCF 17800 X g

Fixed Angle Rotor :

A006-4:

1.5/2.2mL X 12 (16000r/min) RCF: 17800 X g

A008-2: 1.5/2.2mL X 18 (13500r/min) RCF:13550 X g

A011-2: 0.5 mL X 36 (13500r/min) RCF:13250 X g

A011-4: 5mL X 10 (13500r/min) RCF:13250 X g

Power: 100W

Noise <65dB

Power Supply AC110V 60Hz AC220V 50Hz

Dimensions 20X25X30cm